

Detection of nucleic acid cargo in EVs by using Fluorescence Nanoparticle Tracking Analysis

While the characterization of extracellular vesicles (EVs) was mainly focused on proteins and lipids in the past [1-2], it is apparent that interest in the study of nucleic acid cargo in EVs is now increasing. The use of fluorescence nanoparticle tracking analysis (F-NTA) enables the detection of membrane-associated nucleic acids from EVs.



Figure 1.: EV-membrane-associated nucleic acids labeled with $SYBR^{TM}$ Gold.

These experiments used the nucleic acid dye SYBRTM Gold (#S11494, Thermo Fisher Scientific), which has shown a sensitivity of >10x more than ethidium bromide, for the detection of DNA and RNA in denaturing urea, glyoxal and formaldehyde gel [3].

The specificity of the experiment was achieved by a nuclease control aimed at specifically degrading EV-associated nucleic acids after labeling them. However, the use of Benzonase (#70746-3, Novagen) did not allow differentiation between DNA and RNA degradation in this experiment.



Figure 2.: Particle size distribution of SYBR[™] Gold labelled EVs incubated with and without a nuclease prior NTA measurement.

Parameters & Results	Fluorescence λ=488F500nm w/o Nuclease	Fluorescence λ=488F500nm with Nuclease
Detection		
# Positions	11	11
# Cycles	1	1
Camera Parameters		
# Frames	30	30
Sensitivity	95	95
Shutter	200	200
Tracking Parameters		
Tracelength	12	12
Framerate	30 fps	30 fps
Min Area	0	0
Max Area	1000	1000
Min Particle Size	1	1
Max Particle Size	1000	1000
Low Bleach	Yes	Yes
Data Results		
Size (Peak)	82.6	116,5
X50	125.4	140
# Traced Particles	1286	393
Concentration	3x10 ⁷	8x10 ⁶

Table 1: Summary of parameters and results of SYBRTM Gold labeled EVs incubated with and without nuclease.



Application Note



Conclusions:

The reduction of the fluorescent signal after incubation with Benzonase shows that $SYBR^{TM}$ Gold can be used as a specific stain for nucleic acid detection by F-NTA.

It can be assumed that the SYBR[™]-Gold signal remaining after nuclease treatment is due to nucleic acids in the EV lumen that cannot be degraded.

Therefore, SYBR[™]-Gold-activated nucleic acid detection can be used for selective measurement of membrane-associated nucleic acids due to the nuclease treatment.

References:

1. Lötvall et al.:

Minimal experimental requirements for definition of extracellular vesicles and their functions: a position statement from the International Society for Extracellular Vesicles

J Extracell Vesicles 2014; Dec 22:3:26913

2. Théry et al.:

Minimal information for studies of extracellular vesicles 2018 (MISEV2018): a position statement of the International Society for Extracellular Vesicles and update of the MISEV2014 guidelines

J Extracell Vesicles 2018; Nov 23;7(1):1535750

3. Truma et al.:

Characterization of SYBR Gold nucleic acid gel stain: a dye optimized for use with 300-nm ultraviolet transilluminators

Anal Biochem 1999; Mar 15;268(2):278-88

Disclaimer: While we have taken every care to ensure that the information in this document is correct, nothing herein should be construed as implying any representation or warranty as to its accuracy, correctness or completeness of this information and we shall not be held liable for any errors therein or for any damages in connection with the use of this material. Particle Metrix reserves the right to change the content of this material at any time without notice.

Copyright: © May 2024 Particle Metrix. This publication or parts thereof may not be copied or distributed without our express written permission.



Head Office

Particle Metrix GmbH Wildmoos 4 D-82266 Inning / Germany

+49-8143-99172-0 info@particle-metrix.de



US Office

Particle Metrix Inc. Mebane, NC 27302 / USA

+1-919-667-6960 usa@particle-metrix.com



	●

Worldwide Distributors





www.particle-metrix.com